

The JENVIS Dry Eye Report is a new way for eyecare professionals to create a convenient one-page overview of the exam results generated by the Oculus Keratograph 5M, a corneal topographer with a built-in real keratometer and a color camera optimized for external imaging. In addition to the meibography, NIKBUT, tear meniscus height, lipid layer and bulbar redness information from the Keratograph 5M, data collected from other tests, such as tear osmolarity, blink rate and the OSDI Dry Eye Questionnaire can be entered into the report.

The Dry Eye Report, which includes the practice logo, combines screening and consultancy. For the patient, it also provides the doctor's treatment recommendations and easy-to-understand explanations of abbreviations and technical terms. [Click here](#) to see an example of the report.

